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Lykxgf 2-30ml vial compact line for washing, drying, filling and plugging

introduction

The production line is composed of lyqcl series vertical ultrasonic bottle washing machine, lyasmr series tunnel hot air circulation sterilization oven, lykxgf series filling and stoppering machine, lyzgf series capping machine, etc. the production line is simple and beautiful in appearance, small in floor space and convenient in cleaning. It is suitable for the injection product of 2-30ml vial, which can complete the water pouring and super filling of Chinese and Western bottles in drug production. More than 20 process steps such as sonic cleaning, water flushing, air flushing, drying and sterilization, cooling, nitrogen flushing pre-and after filling, plug trimming, plug pressing, cover trimming, cover rolling, etc. The design concept of the compact line is advanced, the structure is reasonable, and the performance is stable. It provides equipment and technical support for the aseptic production of drugs in pharmaceutical plants, and conforms to the relevant requirements of the current national GMP.

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Performance characteristics

- ◆The vertical ultrasonic washing machine is the most popular cleaning equipment in the world, with the best cleaning effect and the best detection effect of visible foreign matters. It has a wide range of application and a low rate of bottle breaking. In the whole cleaning process, the water and gas are all used independent spray needles, each spray needle only fixed with one cleaning medium, a washing and cleaning process is completed to avoid cross contamination and dead angle, which meets the requirements of the new GMP.
- ◆The tunnel type hot air circulation sterilization oven adopts the laminar flow principle and hot air high-speed sterilization process, which can make the container complete the whole production process of preheating, drying and sterilization, cooling to the bottle out under the condition of Class-A laminar flow in the sealed tunnel. It is the most popular drying and sterilization equipment with the best sterilization effect in the world at present, with uniform heat distribution, good pyrogen removal effect, and equipped with mesh belt ultrasonic on-line cleaning, power-off protection and air pressure automatic regulation system.
- ◆The filling and stoppering machine adopts a positive table design structure, and can automatically complete the separation and delivery of the bottle by means of the transfer and delivery of the disk synchronous belt. It is the most popular and the most stable filling and stoppering equipment in the world at present.
- ◆The capping machine adopts the principle of small knife side follow-up rolling capping. The blade can be slightly adjusted in the up and down, and the circumference direction. It is more suitable for domestic packaging materials, and has good adaptability and smooth rolling mouth without damage.
- ◆The whole line adopts advanced intelligent electronic control technology, simple and intuitive human-machine interface, realizes the mechanical and electrical integration design, can be front and back linkage control and single machine operation, ensures the normal operation of the whole line, and can display the machine operation status and fault display, fault self diagnosis, fault analysis and troubleshooting methods, so that the whole production process can achieve automatic control.

设备参数 technical parameters

产品型号 model	LYQCL80+LYASMR620/48+LYKGF8	LYQCL100+LYASMR620/48+LYKGF10	LYQCL120+LYASMR620/55+LYKGF12
适应规格 Adaptive specification	2-30ml国标西林瓶 national standard vials	2-30ml国标西林瓶 national standard vials	2-30ml国标西林瓶 national standard vials
产量 Capacity	16000-12000支/小时 pcs/h	6000-18000支/小时 pcs/h	6000-26000支/小时 pcs/h
澄明度 Clarity	>99%	>99%	>99%
破损率 Breakage rate	≤0.5%(整条线)	≤0.5%(the whole line)	≤0.5%(the whole line)
装量误差 Loading error	±0.5-2%(按药典标准)	±0.5-2%(according to pharmacopoeia standard)	±0.5-2%(according to pharmacopoeia standard)
新鲜水耗量及压力 Fresh water consumption	0.2-0.4Mpa / 0.6-1.0 m ³ /h	0.2-0.4Mpa / 0.6-1.0 m ³ /h	0.2-0.4Mpa / 0.6-1.0 m ³ /h
洁净压缩空气耗量及压力 Clean compressed air consumption and pressure	0.3-0.6Mpa / 35-80 m ³ /h	0.3-0.6Mpa / 35-80 m ³ /h	0.3-0.6Mpa / 35-80 m ³ /h
加塞合格率 And pressure	≥99%	More than 99%	More than 99%
轧盖合格率 Pass rate of capping	≥99%	More than 99%	More than 99%
真空抽速 Vacuum pumping speed	60 m ³ /h (外配)	60 m ³ /h (external distribution)	60 m ³ /h (external distribution)
排风量 Exhaust volume	3500-4800 m ³ /h	3500-4800 m ³ /h	3500-4800 m ³ /h
净重量 Net weight	9000 (kg)	9000(kg)	9000(kg)
总功率 Total power	About 90kw 380V 50Hz	About 93kw 380V 50Hz	About 110KW 380V 50Hz
外形尺寸 Outline size	About 12600 × 2260 × 2445mm	About 12600 × 2260 × 2445mm	About 13200 × 2260 × 2445mm